

Product datasheet for RC204548L1V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

GFAP (NM_002055) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: GFAP (NM_002055) Human Tagged ORF Clone Lentiviral Particle

Symbol: GFAP

Synonyms: ALXDRD

Mammalian Cell None

Selection:

ACCN:

Vector: pLenti-C-Myc-DDK (PS100064)

NM 002055

Tag: Myc-DDK

ORF Size: 1296 bp

ORF Nucleotide

The ORF insert of this clone is exactly the same as(RC204548).

Sequence:

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeg: NM 002055.2

 RefSeq Size:
 3097 bp

 RefSeq ORF:
 1299 bp

 Locus ID:
 2670

 UniProt ID:
 P14136

Cytogenetics: 17q21.31

Domains: filament, filament_head

Protein Families: ES Cell Differentiation/IPS





GFAP (NM_002055) Human Tagged ORF Clone Lentiviral Particle - RC204548L1V

MW: 49.9 kDa

Gene Summary: This gene encodes one of the major intermediate filament proteins of mature astrocytes. It is

used as a marker to distinguish astrocytes from other glial cells during development. Mutations in this gene cause Alexander disease, a rare disorder of astrocytes in the central nervous system. Alternative splicing results in multiple transcript variants encoding distinct

isoforms. [provided by RefSeq, Oct 2008]