

Product datasheet for **RC204548L3V**

GFAP (NM_002055) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	GFAP (NM_002055) Human Tagged ORF Clone Lentiviral Particle
Symbol:	GFAP
Synonyms:	ALXDRD
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_002055
ORF Size:	1296 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204548).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_002055.2
RefSeq Size:	3097 bp
RefSeq ORF:	1299 bp
Locus ID:	2670
UniProt ID:	P14136
Cytogenetics:	17q21.31
Domains:	filament, filament_head
Protein Families:	ES Cell Differentiation/IPS



[View online »](#)

MW: 49.9 kDa

Gene Summary: This gene encodes one of the major intermediate filament proteins of mature astrocytes. It is used as a marker to distinguish astrocytes from other glial cells during development. Mutations in this gene cause Alexander disease, a rare disorder of astrocytes in the central nervous system. Alternative splicing results in multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Oct 2008]