

Product datasheet for **RC204556A1V**

Human MRGPRF (NM_145015) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human MRGPRF (NM_145015) AAV Particle
Tag:	Myc-DDK
Symbol:	MRGPRF
Synonyms:	GPR140; GPR168; MRGF; RTA
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC204556 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCTGGAACTGCTCCTGGGAGGCCCATCCCGCAACAGGAACAAGATGTGCCCTGGCCTGAGCGAGG
CCCCGGAACCTACAGCCGGGGCTTCTGACCATCGAGCAGATCGCGATGCTGCCGCCTCCGGCCGTCAT
GAACTACATCTTCTGCTCCTCTGCCTGTGTGGCCTGGTGGGCAACGGGCTGGTCCTCTGGTTTTTCGGC
TTCTCCATCAAGAGGAACCCCTTCTCCATCTACTTCTGCACCTGGCCAGCGCCGATGTGGGCTACCTCT
TCAGCAAGGCGGTGTTCTCCATCTGAACACGGGGGGCTTCTGGGCACGTTTGCCGACTACATCCGAG
CGTGTGCCGGGTCTGGGGCTCTGCATGTTCTTACCGGCGTGAGCCTCCTGCCGGCCGTCAGCGCCGAG
CGCTGCGCCTCGGTCTCTTCCCGCCTGGTACTGGCGCCGGCGGCCAAGCGCCTGTGCGCCGTGGTGT
GCGCCCTGCTGTGGGTCTGTCCCTCTGGTCACCTGCCTGCACAACCTTCTGCGTGTCTGGGCCG
CGGGGCCCCGGCGCGGCTGCAGGCACATGGACATCTTCTGGGCATCCTCTGTTCTGCTCTGCTGC
CCGCTCATGGTGTGCTGCCCTGCCTGGCCCTCATCTGCACGTGGAGTGCCGGGCCGAGCGGCCAGCGCT
CTGCCAAGCTCAACCAGTCTATCTGGCCATGGTCTCCGTCTTCTGGTGTCTCCATCTACTTAGGGAT
CGACTGGTTCTTCTGGGTCTTCCAGATCCCGGCCCTTCCCGAGTACGTACTGACCTGTGCATC
TGCATCAACAGCAGCGCCAAGCCATCGTCTACTTCTGGCCGGGAGGACAAGTCGACGCGGTGTGGG
AGCCGCTCAGGTGGTCTTCCAGCGGGCCTGCGGGACGCGCTGAGCTGGGGGAGGCCGGGGCAGCAC
GCCCAACACAGTCACCATGGAGATGCAGTGTCCTCCCGGGAACGCCTCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC204556 protein sequence Red=Cloning site Green=Tags(s) MAGNCSWEAHPGNNRNMCPGLSEAPELYSRGFLTIEQIAMLPPPAVMNYIFLLLCLCGLVGNGLVLWFFG FSIKRNPFSIYFLHLASADVGYLFSKAVFSILNTGGFLGTFADYIRSVCRVLGLCMFLTGVSLLPAVSAE RCASVIFPAWYWRRRPKRLSAVVCALLWVLSLLVTCLHNYFCVFLGRGAPGAACRHMDFLGILLFLLCC PLMVLPCALILHVECRARRRQRSAKLNHVILAMVSVFLVSSIYLGIDWFLFWVFQIPAPFPEYVTDLCI CINSSAKPIVYFLAGRDKSQRLWEPLRVVFQRLRDGAELGEAGGSTPNTVTMEMQCPPGNAS TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_145015
ORF Size:	1029 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_145015.2</u>
RefSeq Size:	2308 bp
RefSeq ORF:	1032 bp
Locus ID:	116535
UniProt ID:	<u>Q96AM1</u>
Cytogenetics:	11q13.3
MW:	38.2 kDa