

## Product datasheet for **RC204558L3V**

### CMTM5 (NM\_138460) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	CMTM5 (NM_138460) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CMTM5
Synonyms:	CKLFSF5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_138460
ORF Size:	468 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204558).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_138460.2</a>
RefSeq Size:	1217 bp
RefSeq ORF:	471 bp
Locus ID:	116173
UniProt ID:	<a href="#">Q96DZ9</a>
Cytogenetics:	14q11.2
Protein Families:	Transmembrane
MW:	17.4 kDa



[View online »](#)

**Gene Summary:**

This gene encodes a member of the chemokine-like factor superfamily. This family of genes encodes multi-pass membrane proteins that are similar to both the chemokine and the transmembrane 4 superfamilies of signaling molecules. The encoded protein may exhibit tumor suppressor activity. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2014]