

Product datasheet for **RC204585A1V**

Human TADA1L (TADA1) (NM_053053) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human TADA1L (TADA1) (NM_053053) AAV Particle
Tag:	Myc-DDK
Symbol:	TADA1L
Synonyms:	ADA1; hADA1; HFI1; STAF42; TADA1L
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC204585 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCGACCTTTGTGAGCGAGCTGGAGGCGGCCAAGAAGAACTTAAGCGAGGCCCTGGGGGACAACGTGA
AACAACTACTGGGCTAACCTAAAGCTGTGGTTCAAGCAGAAGATCAGCAAAGAGGAGTTTGACCTTGAAGC
TCATAGACTTCTCACACAGGATAATGTCCATTCTCACAATGATTTCTCCTGGCCATTCTCACGCGTTGT
CAGATTTTGGTTTCTACACCAGATGGTGCTGGATCTTTGCCTTGCCAGGGGGTTCGCAGCAAAACCTG
GAAAACCAAGGAAAGAAAAAGCTTTCTTCTGTTTCGTCAGAAATTTGATCATAGATTCCAGCCTCAAAA
TCCTCTCTCAGGAGCCAGCAATTTGTGGCAAAGGATCCCCAAGATGATGACGACTTGAAACTTTGTGTTCC
CACACAATGATGCTTCCCACTCGAGGCCAGCTTGAAGGAGAATGATAGTGACTGCTTATGAGCATGGGC
TGGACAATGTCACCGAGGAGGCTGTTTCAGCTGTTGTCTATGCTGTGGAGAATCACCTTAAAGATATACT
GACGTCAGTTGTGTCAAGAAGGAAAGCTTATCGGTTACGAGATGGTCATTTTAAATATGCCTTTGGCAGT
AACGTGACCCCGCAGCCATACCTGAAGAATAGTGTAGTAGCTTACAACAATTAATAGAAAGCCCTCCAG
CTTTTACTGCTCCCTGTGCTGGTCAGAATCCAGCTTCTACCCACCCCTGATGATGCTGAGCAGCAGGC
TGCACTCCTGCTGGCATGCTCCGAGAGACTCTACCTGCATCTTGCCTCCGGTGAACATGTACGATCTT
TTTGAAGCTTTGAGGTGCACAGGGAAGTCATCCCTACACATACTGTCTATGCTTTAACATTGAAAGGA
TCATCAGGAACTCTGGCATCCAAATCATGAAGAGCTGCAGCAAGACAAAGTTACCGCCAGCGCTTGGC
AGCCAAGGAGGGGCTTTGCTGTGC

AC**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC204585 protein sequence Red=Cloning site Green=Tags(s) MATFVSELEAAKKNLSEALGDNVKQYWANLKLWFKQKISKEEFDLEAHRLLTQDNVHSHNDFLLAILTRC QILVSTPDGAGSLPWPGGSAAPGKPKGKKLSSVRQKFDHRFPQNPLSGAQQFVAKDPQDDDLKLC HTMMLPTRGQLEGRMIVTAYEHGLDNVTEEAVSAVVYAVENHLKDILTSSVSRRKAYRLRDGHFKYAFGS NVTPQPYLKNSVVAYNNLIESPPAFTAPCAGQNPASHPPDDAEQQAALLLACSGDTLPASLPPVNMVYL FEALQVHREVIPTHTVYALNIERIITKLWHPNHEELQQDKVHRQRLAAKEGLLLC SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_053053
ORF Size:	1005 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_053053.2</u>
RefSeq Size:	2276 bp
RefSeq ORF:	1008 bp
Locus ID:	117143
UniProt ID:	<u>Q96BN2</u>
Cytogenetics:	1q24.1
MW:	37.4 kDa