

Product datasheet for **RC204585L4V**

TADA1L (TADA1) (NM_053053) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TADA1L (TADA1) (NM_053053) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TADA1L
Synonyms:	ADA1; hADA1; HFI1; STAF42; TADA1L
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_053053
ORF Size:	1005 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204585).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053053.2
RefSeq Size:	2276 bp
RefSeq ORF:	1008 bp
Locus ID:	117143
UniProt ID:	Q96BN2
Cytogenetics:	1q24.1
MW:	37.4 kDa


[View online »](#)

Gene Summary:

TADA1L is a protein subunit of the human STAGA complex (SPT3; (MIM 602947)/TAF9 (MIM 600822)/GCN5 (MIM 602301) acetyltransferase complex), which is a chromatin-modifying multiprotein complex (Martinez et al., 2001 [PubMed 11564863]).[supplied by OMIM, Apr 2009]