

Product datasheet for **RC204607A1V**

Human TFIIB (GTF2B) (NM_001514) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human TFIIB (GTF2B) (NM_001514) AAV Particle
Tag:	Myc-DDK
Symbol:	TFIIB
Synonyms:	TF2B; TFIIB
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC204607 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGCGTCTACCAGCCGTTTGGATGCTCTTCCAAGAGTCACATGTCCAAACCATCCAGATGCGATTTTAG
TGGAGGACTACAGAGCCGGTGATATGATCTGCTGAATGTGGCTTGGTTGTAGGTGACCGGGTTATTGA
TGTGGGATCTGAATGGCGAACTTTCAGCAATGACAAAGCAACAAAGATCCATCTCGAGTTGGAGATTCT
CAGAATCCTCTTCTGAGTGATGGAGATTGTCTACCATGATTGGCAAGGGCACAGGAGCTGCAAGTTTTG
ACGAATTTGGCAATTCTAAGTACCAGAATCGGAGAACAATGAGCAGTTCTGATCGGGCAATGATGAATGC
ATTCAAAGAAATCACTACCATGGCAGACAGAATCAATCTACCTCGAAATATAGTTGATCGAACAATAAT
TTATTCAAGCAAGTATATGAACAGAAGAGCCTGAAGGGAAGAGCTAATGATGCTATAGCTTCTGCTTGTC
TCTATATTGCCTGTAGACAAGAAGGGTTCTTAGGACATTTAAAGAAATATGTGCCGTATCACGAATTTCT
TAAGAAAGAAATTGGTCGGTGTTTTAACTTATTTGAAAGCGCTAGAAACCAAGTGTGGATTTGATTACA
ACTGGGGACTTCATGTCCAGGTTCTGTTCCAACCTTTGTCTTCTAAACAAGTACAGATGGCAGCTACAC
ATATAGCCCGTAAAGCTGTGAATTGGACTTGGTTCTGGGAGGAGCCCCATCTGTGGCAGCGGCAGC
TATTTACATGGCCTCACAGGCATCAGCTGAAAAGAGGACCCAAAAAGAAATTGGAGATATTGCTGGTGTT
GCTGATGTTACAATCAGACAGTCCTATAGACTGATCTATCCTCGAGCCCCAGATCTGTTTCTACAGACT
TCAAATTTGACACCCAGTGACAAACTACCACAGCTA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC204607 protein sequence Red=Cloning site Green=Tags(s) MASTRLDALPRVTCPNHPDAILVEDYRAGDMICPECGLVVGDRVIDVGSEWRTFSNDKATKDPSRVGDS QNPLLSGDGLSTMIGKGTGAASFDEFNGSKYQNRRTMSSSDRAMMNAFKEITTMADRINLPNIVDRTNN LFKQVYEQKSLKGRANDAIASACLYIACRQEGVPRTFKEICAVSRISKKEIGRCFKLILKALETSDLIT TGDMSRFCSNLCLPKQVQMAATHIARKAVELDLVPGRSPISVAAAAIYMASQASAEKRTQKEIGDIAGV ADVTIRQSYRLIYPAPDLFPTDFKFDTPVDKLPQL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001514
ORF Size:	948 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001514.3</u>
RefSeq Size:	1651 bp
RefSeq ORF:	951 bp
Locus ID:	2959
UniProt ID:	<u>Q00403</u>
Cytogenetics:	1p22.2
MW:	34.8 kDa