

Product datasheet for **RC204610A1V**

Human TMEM156 (NM_024943) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human TMEM156 (NM_024943) AAV Particle
Tag:	Myc-DDK
Symbol:	TMEM156
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC204610 representing NM_024943 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGACAAAAACAGCCCTCCTTAAATTATTTGTGGCAATAGTGATCACATTCATTTTAATTTTGCCGGAAT
ATTTCAAGACACCGAAAGAAAGAACATTGGAGCTATCATGTCTGGAAGTGTGTTTGCAATCTAATTTTAC
CTATTCACCTCCTCCTTAAATTTTCTTTGTGACTTTCTGCAACCAGTAAGGGAAACTCAGATTATC
ATGAGAATCTTTCTAAATCCCTCAATTTTCGTAACCTCACCAGGACTTGCCAAGACATCACAGGTGAAT
TAAAAATGTGCTCCTCGTGTGTTGTTGTGAGCCTAAAGGAAACATGGATTTTATTTCTCAGGAACAAAC
ATCAAAAGTTCTTATCAGGAGAGGATCAATGGAAGTGAAAGCAAATGATTTTCATTCACCTTGTCAGCAC
TTTAACTTCAGTGTAGCTCCTCTGGTTGACCACTTGGAAGGAATATAACACTACCTGTCATCTAAAAAACC
ACACTGGAAGATCAACAATCATGGAGGATGAGCCAAGCAAGGAGAAATCGATAAACTACACTTGTAGAAT
CATGGAATACCGAATGATTGTATACACATTTCTTTGCACCTGGAGATGGATATAAAAAATATCACTTGT
TCCATGAAGATCACTTGGTATATTTTAGTTCTATTAGTTTTATATTTTGTATCATCCTCACTATCCGCA
AAATACTTGAAGGCCAGAGAAGAGTGCAAAAGTGGCAGAGTCATAGAGACAAACCTACATCTGTTCTCTT
AAGAGGAAGTGATTCGGAGAACTGAGAGCATTGAATGTGCAGGTTCTTTCAGCAGAGACCACGCAGAGG
CTGCCTTTGGATCAAGTCCAGGAAGTGCTTCCCCCAATTCAGAACTA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC204610 representing NM_024943 Red=Cloning site Green=Tags(s) MTKTALLKLFVAIVITFILILPEYFKTPKERTLELSCLEVCLQSNFTYSLSSLNFSFVTFLQPVRETQII MRIFLNPSNFRNFRTRCQDITGEFKMCSSCLVCEPKGNDMDFISQEOTSKVLIRRGSMVEKANDFHSPCQH FNFSVAPLVDHLEEYNTTCHLKNHTGRSTIMEDEPSKEKSINYTCRIMEYPNDCIHISLHLEMDIKNITC SMKITWYILVLLVFIFLIILTIRKILEGQRRVQKWQSHRDKPTSVLLRGSDSEKLRALNVQVLSAETTQR LPLDQVQEVLPPIPEL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_024943
ORF Size:	888 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_024943.1</u>
RefSeq Size:	1880 bp
RefSeq ORF:	891 bp
Locus ID:	80008
UniProt ID:	<u>Q8N614</u>
Cytogenetics:	4p14
MW:	34.1 kDa