

Product datasheet for RC204610L4V

TMEM156 (NM_024943) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TMEM156 (NM_024943) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TMEM156
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_024943
ORF Size:	888 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204610).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024943.1
RefSeq Size:	1880 bp
RefSeq ORF:	891 bp
Locus ID:	80008
UniProt ID:	Q8N614
Cytogenetics:	4p14
Protein Families:	Transmembrane
MW:	34.1 kDa


[View online »](#)