

Product datasheet for **RC204629L4V**

Alpha 1 Acid Glycoprotein (ORM1) (NM_000607) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Alpha 1 Acid Glycoprotein (ORM1) (NM_000607) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Alpha 1 Acid Glycoprotein
Synonyms:	AGP-A; AGP1; HEL-S-153w; ORM
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_000607
ORF Size:	603 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204629).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_000607.1
RefSeq Size:	847 bp
RefSeq ORF:	606 bp
Locus ID:	5004
UniProt ID:	P02763
Cytogenetics:	9q32
Domains:	lipocalin
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)

MW: 23.5 kDa

Gene Summary: This gene encodes a key acute phase plasma protein. Because of its increase due to acute inflammation, this protein is classified as an acute-phase reactant. The specific function of this protein has not yet been determined; however, it may be involved in aspects of immunosuppression. [provided by RefSeq, Jul 2008]