

Product datasheet for **RC204650L1V**

Nuclear Factor 1 (NFIA) (NM_005595) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nuclear Factor 1 (NFIA) (NM_005595) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Nuclear Factor 1
Synonyms:	BRMUTD; CTF; NF-I/A; NF1-A; NFI-A; NFI-L
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_005595
ORF Size:	1494 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204650).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_005595.1
RefSeq Size:	9403 bp
RefSeq ORF:	1497 bp
Locus ID:	4774
UniProt ID:	Q12857
Cytogenetics:	1p31.3
Protein Families:	Transcription Factors
MW:	54.6 kDa



[View online »](#)

Gene Summary:

This gene encodes a member of the NF1 (nuclear factor 1) family of transcription factors. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2011]