

Product datasheet for **RC204670L4V**

SLC35A1 (NM_006416) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SLC35A1 (NM_006416) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLC35A1
Synonyms:	CDG2F; CMPST; CST; hcST
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_006416
ORF Size:	1011 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204670).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006416.2
RefSeq Size:	1918 bp
RefSeq ORF:	1014 bp
Locus ID:	10559
UniProt ID:	P78382
Cytogenetics:	6q15
Domains:	Nuc_sug_transp
Protein Families:	Transmembrane



[View online »](#)

MW: 36.8 kDa

Gene Summary: The protein encoded by this gene is found in the membrane of the Golgi apparatus, where it transports nucleotide sugars into the Golgi. One such nucleotide sugar is CMP-sialic acid, which is imported into the Golgi by the encoded protein and subsequently glycosylated. Defects in this gene are a cause of congenital disorder of glycosylation type 2F (CDG2F). Two transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Dec 2009]