

Product datasheet for **RC204691L3V**

Macrophage Inflammatory Protein 3 alpha (CCL20) (NM_004591) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Macrophage Inflammatory Protein 3 alpha (CCL20) (NM_004591) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Macrophage Inflammatory Protein 3 alpha
Synonyms:	CKb4; Exodus; LARC; MIP-3-alpha; MIP-3a; MIP3A; SCYA20; ST38
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_004591
ORF Size:	285 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204691).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_004591.1
RefSeq Size:	851 bp
RefSeq ORF:	291 bp
Locus ID:	6364
UniProt ID:	P78556
Cytogenetics:	2q36.3
Protein Families:	Druggable Genome, Secreted Protein



[View online »](#)

Protein Pathways: Chemokine signaling pathway, Cytokine-cytokine receptor interaction

MW: 10.7 kDa

Gene Summary: This antimicrobial gene belongs to the subfamily of small cytokine CC genes. Cytokines are a family of secreted proteins involved in immunoregulatory and inflammatory processes. The CC cytokines are proteins characterized by two adjacent cysteines. The protein encoded by this gene displays chemotactic activity for lymphocytes and can repress proliferation of myeloid progenitors. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2014]