

Product datasheet for **RC204696L2V**

CASQ1 (NM_001231) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CASQ1 (NM_001231) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CASQ1
Synonyms:	CASQ; PDIB1; VMCQA
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_001231
ORF Size:	1170 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204696).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001231.3 , NP_001222.2
RefSeq Size:	1993 bp
RefSeq ORF:	1191 bp
Locus ID:	844
UniProt ID:	P31415
Cytogenetics:	1q23.2
Domains:	Calsequestrin
MW:	44.5 kDa


[View online »](#)

Gene Summary:

This gene encodes the skeletal muscle specific member of the calsequestrin protein family. Calsequestrin functions as a luminal sarcoplasmic reticulum calcium sensor in both cardiac and skeletal muscle cells. This protein, also known as calmitine, functions as a calcium regulator in the mitochondria of skeletal muscle. This protein is absent in patients with Duchenne and Becker types of muscular dystrophy. [provided by RefSeq, Jun 2013]