

Product datasheet for **RC204702A1V**

Human Pregnancy specific beta 1 glycoprotein 11 (PSG11) (NM_203287) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Human Pregnancy specific beta 1 glycoprotein 11 (PSG11) (NM_203287) AAV Particle
Tag: Myc-DDK
Symbol: Pregnancy specific beta 1 glycoprotein 11
Synonyms: PSBG-11; PSBG-13; PSG13; PSG14
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >RC204702 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCACGCAGCAGAGATCATGGGGCCCCTCTCAGCCCCTCCCTGCACAGAGCACATCAAATGGAAGGGGC
 TCCTGCTCACAGTGGAGACTCCCAAGCCCTCCATCTCCAGCAGCAACTTAAACCCAGGGAGGCCATGGA
 GACTGTGATCTTAACCTGTAATCCTGAGACTCCGGACGCAAGCTACCTGTGGTGGATGAATGGTCAGAGC
 CTCCCTATGACTCATAGGATGCAGCTGTCTGAAACCAACAGGACCCTCTTTCTATTTGGTGCACAAAGT
 ATACTGCAGGACCCTATGAATGTGAAATATGGAACCTCAGGGAGTGCCAGCCGAGTGACCCAGTCCACC
 GAATCTCTCCATGGTCCAGACCTCCCCAGAATTTTCCCTTCAGTCACCTTCTACTATTCAGGAGAGAAC
 CTCGACTTGTCTGCTTCGCAGACTCTAACCCACCAGCACAGTATTCTTGGACAATTAATGGGAAGTTTC
 AGCTATCAGGACAAAAGCTCTTTATCCCTCAAATTACTCAAAGCATAATGGGCTCTATGCTTGCTCTGC
 TCGTAACCTCAGCCACTGGCGAGGAAAGCTCCACATCCTTGACAATCAGAGTCATTGCTCTCCAGGATTA
 GGAACTTTTGCTTTCAATAATCCAACG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC204702 protein sequence
 Red=Cloning site Green=Tags(s)

MHAAEIMGPLSAPPCTEHIKWKGLLLTVETPKPSISSNLPREAMETVILTCNPETPDASYLWWMNGQS
 LPMTHRMQLSETNRTLFLFGVTKYTAGPYECEIWNSSASRSDPVTLLNLLHGPDLPRIFPSVTSYSGEN
 LDLSCFADSNPPAQYSWTINGKFKLSGQKLFIPQITPKHNGLYACARSATGEESSLSLTIRVIAPPGL
 GTFAFNNPT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)

Species:	Human
Serotype:	AAV-2
ACCN:	NM_203287
ORF Size:	660 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_203287.1
RefSeq Size:	1243 bp
RefSeq ORF:	642 bp
Locus ID:	5680
UniProt ID:	Q9UQ72
Cytogenetics:	19q13.31
MW:	24 kDa