

Product datasheet for **RC204713L3V**

Apolipoprotein F (APOF) (NM_001638) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Apolipoprotein F (APOF) (NM_001638) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Apolipoprotein F
Synonyms:	Apo-F; LTIP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001638
ORF Size:	978 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204713).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001638.2
RefSeq Size:	1720 bp
RefSeq ORF:	981 bp
Locus ID:	319
UniProt ID:	Q13790
Cytogenetics:	12q13.3
Protein Families:	Druggable Genome, Secreted Protein
MW:	35.4 kDa



[View online »](#)

Gene Summary:

The product of this gene is one of the minor apolipoproteins found in plasma. This protein forms complexes with lipoproteins and may be involved in transport and/or esterification of cholesterol. [provided by RefSeq, Jul 2008]