

Product datasheet for RC204727L2

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com

CN: techsupport@origene.cn

OriGene Technologies, Inc.

CLIC2 (NM_001289) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: CLIC2 (NM_001289) Human Tagged Lenti ORF Clone

Tag: mGFP Symbol: CLIC2

Synonyms: CLCNL2; CLIC2b; MRXS32; XAP121

Mammalian Cell None

Selection:

Vector:pLenti-C-mGFP (PS100071)E. coli Selection:Chloramphenicol (34 ug/mL)


ORF Nucleotide The ORF insert of this clone is exactly the same as(RC204727).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





^{*} The last codon before the Stop codon of the ORF.

ACCN: NM_001289

ORF Size: 741 bp



CLIC2 (NM_001289) Human Tagged Lenti ORF Clone - RC204727L2

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 001289.4</u>

 RefSeq Size:
 2694 bp

 RefSeq ORF:
 744 bp

 Locus ID:
 1193

 UniProt ID:
 015247

Cytogenetics: Xq28

Protein Families: Druggable Genome, Ion Channels: Other

MW: 28.4 kDa


Gene Summary: This gene encodes a chloride intracellular channel protein. Chloride channels are a diverse

group of proteins that regulate fundamental cellular processes including stabilization of cell membrane potential, transepithelial transport, maintenance of intracellular pH, and regulation of cell volume. This protein plays a role in inhibiting the function of ryanodine receptor 2. A mutation in this gene is the cause of an X-linked form of cognitive disability.

[provided by RefSeg, Jul 2017]



Product images:



Circular map for RC204727L2