

Product datasheet for **RC204734A1V**

Human LTBR (NM_002342) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human LTBR (NM_002342) AAV Particle
Tag:	Myc-DDK
Symbol:	LTBR
Synonyms:	D12S370; LT-BETA-R; TNF-R-III; TNFCR; TNFR-RP; TNFR2-RP; TNFR3; TNFRSF3
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC204734 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGCTCCTGCCTTGGGCCACCTCTGCCCCGGCCTGGCCTGGGGCCTCTGGTGCTGGGCCTCTTCGGGC
TCCTGGCAGCATCGCAGCCCCAGGCGGTGCCTCCATATGCGTCGGAGAACCAGACCTGCAGGGACCAGGA
AAAGGAATACTATGAGCCCCAGCACCGCATCTGCTGCTCCCGCTGCCCGCCAGGCACCTATGTCTCAGCT
AAATGTAGCCGCATCCGGGACACAGTTTGTGCCACATGTGCCGAGAATTCCTACAACGAGCACTGGAAC
ACCTGACCATCTGCCAGCTGTGCCGCCCTGTGACCCAGTGATGGCCTCGAGGAGATTGCCCCCTGCAC
AAGCAAACGGAAGACCCAGTGCCGCTGCCAGCCGGAATGTTCTGTGCTGCCTGGGCCCTCGAGTGTACA
CACTGCCGAGCTACTTTCTGACTGCCCGCTGGCACTGAAGCCGAGCTCAAAGATGAAGTTGGGAAGGGTA
ACAACCACTGCGTCCCCTGCAAGGCCGGGCACTTCCAGAATACCTCCTCCCCAGCGCCCGCTGCCAGCC
CCACACCAGGTGTGAGAACCAAGGTCTGGTGGAGGCAGCTCCAGGCACTGCCAGTCCGACACAACCTGC
AAAAATCCATTAGAGCCACTGCCCCAGAGATGTCAGGAACCATGCTGATGCTGGCCGTTCTGCTGCCAC
TGGCCTTCTTTCTGCTCCTTGCCACCGTCTTCTCCTGCATCTGGAAGAGCCACCTTCTCTCTGCAGGAA
ACTGGGATCGTGCTCAAGAGGCGTCCGAGGGAGAGGGACCAATCCTGTAGCTGGAAGCTGGGAGCCT
CCGAAGGCCCATCCATACTTCCCTGACTTGGTACAGCCACTGCTACCCATTTCTGGAGATGTTTCCCGAG
TATCCACTGGGCTCCCCGAGCCCCAGTTTGGAGGCAGGGGTGCCGAACAGCAGAGTCTCTGGACCT
GACCAGGGAGCCGAGTTGGAACCCGGGAGCAGAGCCAGGTGGCCACGGTACCAATGGCATTATGTC
ACCGGCGGGTCTATGACTATCACTGGCAACATCTACATCTACAATGGACCAGTACTGGGGGGACCACCGG
GTCTGGAGACCTCCAGCTACCCCCGAACCTCCATACCCATTCCCAGAGGGGGACCTGGCCCTCC
CGGGCTCTCTACACCCACCAGGAAGATGGCAAGGCTTGGCACCTAGCGGAGACAGAGCACTGTGGTGCC
ACACCCCTAACAGGGGCCCAAGGAACCAATTTATACCCATGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC204734 protein sequence Red=Cloning site Green=Tags(s) MLLPWATSAPGLAWGPLVLGLFGLLAASQPQAVPPYASENQTCRDQEKEYYEPQHRICCSRCPPGTYYVSA KCSRIRDTCATCAENSYNEHWNLTICQLCRPCDPVMGLEEIPCTSKRKTQCRCQPGMFCAAWALECT HCELLSDCPPGTEAEKDEVGKGNHCVCKAGHFQNTSSPSARCQPHTRCENQGLVEAAPGTAQSDTTC KNPLEPLPPEMSGTMLMLAVLLPLAFFLLLATVFSIWKSHPSLCRKLGSLKRRPQGEQGNPVAGSWEPE PKAHPYFDDLQPLLPISGDVSPVSTGLPAAPVLEAGVPQQQSPLDLTREPQLEPGEQSQVAHGTINGIHV TGGSMITIGNIYINGPVLGGPPGPDLPATPEPPYPIPEEGDPGPPGLSTPHQEDGKAWHLAETEHCGA TPSNRGRPNQFITHD TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_002342
ORF Size:	1305 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_002342.1</u>
RefSeq Size:	2283 bp
RefSeq ORF:	1308 bp
Locus ID:	4055
UniProt ID:	<u>P36941</u>
Cytogenetics:	12p13.31
MW:	46.7 kDa