

Protein Sequence:	>RC204746 protein sequence Red=Cloning site Green=Tags(s) MLLFLLSALVLLTQPLGYLEAGMKTYSHRTMPSACTLVMCSSVESGLPGRDGRDREGPRGEKGDPLPG AAGQAGMPGQAGPVGPKGDNDSVGEPPGPKGDTGPSGPPGPPGVPAGREGPLGKQGNIGPQKPGPKGE AGPKGEVGAPGMQGSAGARGLAGPKGERGVPGERGVPNGTGAAGSAGAMGPQGSFGARGPPGLKGDGIP GDKGAKGESGLPDVASLRQQVEALQGQVQHLQAAFSQYKKVELFPNGQSVGEKIFKTAGFVKPFTEAQLL CTQAGGQLASPRSAANAALQQLVYAKNEAFLSMTDSKTEGKFTYPTGESLVYSNWAPGKPNDDGGSED CVEIFTNGKWNDRACGEKRLVVECF TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_003019
ORF Size:	1125 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_003019.4
RefSeq Size:	1299 bp
RefSeq ORF:	1128 bp
Locus ID:	6441
UniProt ID:	P35247
Cytogenetics:	10q22.3
MW:	37.7 kDa