

Product datasheet for **RC204773L4V**

CDA08 (ITFG1) (NM_030790) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	CDA08 (ITFG1) (NM_030790) Human Tagged ORF Clone Lentiviral Particle
Symbol:	CDA08
Synonyms:	2310047C21Rik; CDA08; LNKN-1; TIP
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_030790
ORF Size:	1836 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204773).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_030790.2
RefSeq Size:	2261 bp
RefSeq ORF:	1839 bp
Locus ID:	81533
UniProt ID:	Q8TB96
Cytogenetics:	16q12.1
Domains:	FG-GAP
Protein Families:	Secreted Protein, Transmembrane



[View online »](#)

MW: 68.1 kDa

Gene Summary: Modulator of T-cell function. Has a protective effect in graft versus host disease model (By similarity).[UniProtKB/Swiss-Prot Function]