

Product datasheet for **RC204782L4V**

ARHGAP15 (NM_018460) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	ARHGAP15 (NM_018460) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ARHGAP15
Synonyms:	BM046
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_018460
ORF Size:	1425 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204782).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_018460.2
RefSeq Size:	1774 bp
RefSeq ORF:	1428 bp
Locus ID:	55843
UniProt ID:	Q53QZ3
Cytogenetics:	2q22.2-q22.3
MW:	54.5 kDa



[View online »](#)

Gene Summary:

RHO GTPases (see ARHA; MIM 165390) regulate diverse biologic processes, and their activity is regulated by RHO GTPase-activating proteins (GAPs), such as ARHGAP15 (Seoh et al., 2003 [PubMed 12650940]).[supplied by OMIM, Mar 2008]