

Product datasheet for **RC204792A1V**

Human NGDN (NM_001042635) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human NGDN (NM_001042635) AAV Particle
Tag:	Myc-DDK
Symbol:	NGDN
Synonyms:	C14orf120; CANu1; LCP5; Ipd-2; NGD
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC204792 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCGGCGCTGGGGGTGCTGGAGTCCGACCTGCCAAGTGCCGTGACACTTCTGAAAAATCTCCAGGAGC
AAGTGATGGCTGTAAGTGCACAAGTGAATCACTGACACAAAAGTTCAAGCTGGTGCCTATCCTACAGA
AAAGGGTCTCAGCTTCTTGAAGTGAAGACCAGCTGCTGCTCATGTACCTTATGGATTTGACCCACCTC
ATTCTGGACAAAGCCTCAGGAGGATCTCTTCAGGGACATGATGCAGTTTTGAGACTGGTAGAGATTCGAA
CGGTTTTGAAAAGCTTCGTCCCTTGGACCAAAAGCTGAAGTATCAAATTGACAAGCTGATCAAGACTGC
AGTGACAGGCAGCCTTAGTGAGAATGACCCACTTCGTTTTTAAGCCTCATCCAGCAATATGATGAGCAAG
TTGAGCTCTGAGGATGAGGAGGAAGATGAAGCAGAAGATGACCACTGAGGCTTCAGGGAAGAAATCTG
TGAAGGGAGTGTCTAAGAAATATGTTCTCCACGCTTGGTTCCAGTACATTATGATGAAACAGAAGCTGA
GCGGGAGAAGAAGCGTCTAGAACGAGCCAAGAGACGGGCATTGAGCAGCTCTGTCATTCTGAACTTAAG
GAGCAGTACTCAGATGCTCCAGAGGAAATCCGTGATGCTCGGCATCCCATGTTACCCGCCAGAGTCAGG
AGGACCAACACAGGATTAAGTATGAGGAGAGCATGATGGTGCGTTTGAGCGTCAGTAAGCGAGAGAAAGG
ACGGCGAAAACGAGCAAATGTCATGAGCTCACAACCTTATTCCCTTACACACTTCAGTGACATCAGTGCT
TTGACAGGGGAACTGTTTCATCTTGATGAGGATCAGAATCCTATTAAGAAGCGGAAGAAGATACCTCAGA
AAGTCCGAAGAAAAAGTTTTCGGAGCGGCCG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC204792 protein sequence Red=Cloning site Green=Tags(s) MAALGVLESDLPASVTLLKNLQEQVMAVTAQVKSLTQKVQAGAYPTEKGLSFLEVKDQLLLMYLMDLTHL ILDKASGGSLQGHDAVLRLVEIRTVLEKLRPLDQKLKYQIDKL IKTAVTGSLSENDPLRFKPHPSNMMSK LSSEDEEEDEAEDDQSEASGKKS VKGVSKKYVPPRLVPVHYDETEAEREKKRLERAKRRALSSSVIRELK EQYSDAPEEIRDARHPHVTRQSQEDQHRINYEESMMVRLSVSKREKGRRKRANVMSSQLHSLTHFSDISA LTGGTVHLDEDQNPIKKRKKIPQKGRKKKGFRRRR TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001042635
ORF Size:	945 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_001042635.1, NP_001036100.1</u>
RefSeq Size:	1142 bp
RefSeq ORF:	948 bp
Locus ID:	25983
UniProt ID:	<u>Q8NEJ9</u>
Cytogenetics:	14q11.2
MW:	35.9 kDa