

## Product datasheet for **RC204793L1V**

### PLAC1 (NM\_021796) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	PLAC1 (NM_021796) Human Tagged ORF Clone Lentiviral Particle
Symbol:	PLAC1
Synonyms:	CT92; OOSP2B; OOSP2L
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_021796
ORF Size:	636 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204793).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_021796.3</a>
RefSeq Size:	1126 bp
RefSeq ORF:	639 bp
Locus ID:	10761
UniProt ID:	<a href="#">Q9HBJ0</a>
Cytogenetics:	Xq26.3
Protein Families:	Secreted Protein
MW:	23.4 kDa



[View online »](#)

**Gene Summary:** May play a role in placental development.[UniProtKB/Swiss-Prot Function]