

Product datasheet for RC204803L2

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Osteopontin (SPP1) (NM_001040058) Human Tagged Lenti ORF Clone

Product data:

Product Type: Expression Plasmids

Product Name: Osteopontin (SPP1) (NM_001040058) Human Tagged Lenti ORF Clone

Tag: mGFP

Symbol: Osteopontin

Synonyms: BNSP; BSPI; ETA-1; OPN

Mammalian Cell None

Selection:

Vector: pLenti-C-mGFP (PS100071)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide


The ORF insert of this clone is exactly the same as(RC204803).

Sequence:

Restriction Sites: Sgfl-Mlul

Cloning Scheme:





st The last codon before the Stop codon of the ORF.

ACCN: NM_001040058

ORF Size: 942 bp





Osteopontin (SPP1) (NM_001040058) Human Tagged Lenti ORF Clone - RC204803L2

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method: 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

RefSeq: <u>NM 001040058.1</u>

 RefSeq Size:
 1641 bp

 RefSeq ORF:
 945 bp

 Locus ID:
 6696

 UniProt ID:
 P10451

Cytogenetics: 4q22.1

Protein Families: Druggable Genome, Secreted Protein

Protein Pathways: ECM-receptor interaction, Focal adhesion, Toll-like receptor signaling pathway

MW: 35.4 kDa

Gene Summary: The protein encoded by this gene is involved in the attachment of osteoclasts to the


mineralized bone matrix. The encoded protein is secreted and binds hydroxyapatite with high

affinity. The osteoclast vitronectin receptor is found in the cell membrane and may be involved in the binding to this protein. This protein is also a cytokine that upregulates expression of interferon-gamma and interleukin-12. Several transcript variants encoding

different isoforms have been found for this gene. [provided by RefSeq, Oct 2011]



Product images:



Circular map for RC204803L2