

Product datasheet for **RC204809L1V**

TSPAN2 (NM_005725) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	TSPAN2 (NM_005725) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TSPAN2
Synonyms:	NET3; TSN2; TSPAN-2
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_005725
ORF Size:	663 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204809).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_005725.3
RefSeq Size:	3213 bp
RefSeq ORF:	666 bp
Locus ID:	10100
UniProt ID:	O60636
Cytogenetics:	1p13.2
Domains:	transmembrane4
Protein Families:	Transmembrane



[View online »](#)

MW: 24.1 kDa

Gene Summary: The protein encoded by this gene is a member of the transmembrane 4 superfamily, also known as the tetraspanin family. Most of these members are cell-surface proteins that are characterized by the presence of four hydrophobic domains. The proteins mediate signal transduction events that play a role in the regulation of cell development, activation, growth and motility. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Apr 2015]