

## Product datasheet for **RC204825L2V**

### NMNAT1 (NM\_022787) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	NMNAT1 (NM_022787) Human Tagged ORF Clone Lentiviral Particle
Symbol:	NMNAT1
Synonyms:	LCA9; NMNAT; PNAT1; SHILCA
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_022787
ORF Size:	837 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204825).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_022787.2</a>
RefSeq Size:	3781 bp
RefSeq ORF:	840 bp
Locus ID:	64802
UniProt ID:	<a href="#">Q9HAN9</a>
Cytogenetics:	1p36.22
Domains:	CTP_transf_2
Protein Pathways:	Metabolic pathways, Nicotinate and nicotinamide metabolism



[View online »](#)

**MW:** 31.9 kDa

**Gene Summary:** This gene encodes an enzyme which catalyzes a key step in the biosynthesis of nicotinamide adenine dinucleotide (NAD). The encoded enzyme is one of several nicotinamide nucleotide adenylyltransferases, and is specifically localized to the cell nucleus. Activity of this protein leads to the activation of a nuclear deacetylase that functions in the protection of damaged neurons. Mutations in this gene have been associated with Leber congenital amaurosis 9. Alternative splicing results in multiple transcript variants. Pseudogenes of this gene are located on chromosomes 1, 3, 4, 14, and 15. [provided by RefSeq, Jul 2014]