

## Product datasheet for **RC204832L4V**

### **TOMM40L (NM\_032174) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	TOMM40L (NM_032174) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TOMM40L
Synonyms:	TOMM40B
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_032174
ORF Size:	924 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204832).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_032174.4</a>
RefSeq Size:	2927 bp
RefSeq ORF:	927 bp
Locus ID:	84134
UniProt ID:	<a href="#">Q969M1</a>
Cytogenetics:	1q23.3
Domains:	Euk_porin
Protein Families:	Druggable Genome, Ion Channels: Other



[View online »](#)

**Protein Pathways:** Amyotrophic lateral sclerosis (ALS)

**MW:** 33.9 kDa

**Gene Summary:** Potential channel-forming protein implicated in import of protein precursors into mitochondria.[UniProtKB/Swiss-Prot Function]