

Product datasheet for **RC204852A1V**

Human Granzyme A (GZMA) (NM_006144) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human Granzyme A (GZMA) (NM_006144) AAV Particle
Tag:	Myc-DDK
Symbol:	Granzyme A
Synonyms:	CTLA3; HFSP
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC204852 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGAGGAACCTCTATAGATTTCTGGCATCCTCTCTCAGTTGTCGTTTCTCTCTGCTAATTCCTGAAG
ATGCTGTGAAAAAATTATTGGAGGAAATGAAGTAACCTCATTCAAGACCTACATGGTCCTACTTAG
TCTTGACAGAAAAACCATCTGTGCTGGGGCTTTGATTGCAAAAGACTGGGTGTTGACTGCAGCTCACTGT
AACTTGAACAAAAGGTCCAGGTCATTCTTGGGGCTCACTCAATAACCAGGGAAGAGCCAACAAAACAGA
TAATGCTTGTTAAGAAAGAGTTCCCTATCCATGCTATGACCCAGCCACACGCGAAGGTGACCTTAACT
TTTACAGCTGACGAAAAAGCAAAAATTAACAAATATGTGACTATCCTTCATCTACCTAAAAAGGGGGAT
GATGTGAAACCAGGAACCATGTGCCAAGTTGCAGGGTGGGGCAGGACTCACAATAGTGCATCTTGGTCCG
ATACTCTGAGAGAAGTCAATATCACCATCATAGACAGAAAAGTCTGCAATGATCGAAATCACTATAATTT
TAACCCTGTGATTGGAATGAATATGGTTTGTGCTGGAAGCCTCCGAGGTGGAAGAGACTCGTGAATGGA
GATTCTGGAAGCCCTTTGTTGTGCGAGGGTGTTCGAGGGGTCACTTCCTTTGGCCTTGAAAAATAAT
GCGGAGACCCTCGTGGGCTGGTGTCTATTTCTCTCAAAGAAACACCTCACTGGATAATTATGAC
TATCAAGGGAGCAGTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC204852 protein sequence Red=Cloning site Green=Tags(s) MRNSYRFLASSLSVVVSLLLIPEDVCEKIIGGNEVTPHSRPMVLLSLDRKTICAGALIAKDWLTAHC NLNKRSQLGAHSITREEPTKQIMLVKKEFPYPCYDPATREGDLKLLQLTEKAKINKYVTILHLPKKG DVKPGTMCQVAGWGRTHNSASWSDTLREVNITIIDRKVCNDRNHYNFNPVIGMNMVCAGSLRGG RDSNG DSGSPLLCEGVFRGVTSFGLENKCGDPRGPGVYILLSKKHLNWIIMTIKGA TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_006144
ORF Size:	786 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	<u>NM_006144.2</u>
RefSeq Size:	913 bp
RefSeq ORF:	789 bp
Locus ID:	3001
UniProt ID:	<u>P12544</u>
Cytogenetics:	5q11.2
MW:	29 kDa