

## Product datasheet for **RC204859L1V**

### SCY1 like 3 (SCYL3) (NM\_020423) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SCY1 like 3 (SCYL3) (NM_020423) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SCY1 like 3
Synonyms:	PACE-1; PACE1
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_020423
ORF Size:	2064 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204859).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_020423.4</a> , <a href="#">NP_065156.4</a>
RefSeq Size:	2928 bp
RefSeq ORF:	2067 bp
Locus ID:	57147
UniProt ID:	<a href="#">Q8IZE3</a>
Cytogenetics:	1q24.2
Domains:	S_TKc
Protein Families:	Druggable Genome, Protein Kinase


[View online »](#)

**MW:** 76.7 kDa

**Gene Summary:** This gene encodes a protein with a kinase domain and four HEAT repeats. The encoded protein interacts with the C-terminal domain of ezrin, an ERM protein, and may play a role in cell adhesion and migration. Alternative splicing results in multiple transcript variants encoding multiple isoforms. [provided by RefSeq, Jun 2012]