

Product datasheet for **RC204865L3V**

SLC5A6 (NM_021095) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	SLC5A6 (NM_021095) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLC5A6
Synonyms:	NERIB; SMVT
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_021095
ORF Size:	1905 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204865).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_021095.1
RefSeq Size:	3272 bp
RefSeq ORF:	1908 bp
Locus ID:	8884
UniProt ID:	Q9Y289
Cytogenetics:	2p23.3
Domains:	SSF
Protein Families:	Transmembrane



[View online »](#)

MW: 68.6 kDa

Gene Summary: Transports pantothenate, biotin and lipoate in the presence of sodium.[UniProtKB/Swiss-Prot Function]