

Product datasheet for **RC204875L4V**

PPP2R5C (NM_178588) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | PPP2R5C (NM_178588) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | PPP2R5C |
| Synonyms: | B56G; MGC23064; PR61G |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_178588 |
| ORF Size: | 1152 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC204875). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_178588.1 , NP_848703.1 |
| RefSeq Size: | 1460 bp |
| RefSeq ORF: | 1154 bp |
| Locus ID: | 5527 |
| Cytogenetics: | 14q32.31 |
| Protein Families: | Druggable Genome, Phosphatase |
| Protein Pathways: | Oocyte meiosis, Wnt signaling pathway |
| MW: | 44.9 kDa |



[View online »](#)

Gene Summary:

The product of this gene belongs to the phosphatase 2A regulatory subunit B family. Protein phosphatase 2A is one of the four major Ser/Thr phosphatases, and it is implicated in the negative control of cell growth and division. It consists of a common heteromeric core enzyme, which is composed of a catalytic subunit and a constant regulatory subunit, that associates with a variety of regulatory subunits. The B regulatory subunit might modulate substrate selectivity and catalytic activity. This gene encodes a gamma isoform of the regulatory subunit B56 subfamily. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008]