

# Product datasheet for RC204877L3

### OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

# GRO beta (CXCL2) (NM\_002089) Human Tagged Lenti ORF Clone

### **Product data:**

**Product Type:** Expression Plasmids

Product Name: GRO beta (CXCL2) (NM\_002089) Human Tagged Lenti ORF Clone

Tag: Myc-DDK
Symbol: GRO beta

Synonyms: CINC-2a; GRO2; GROb; MGSA-b; MIP-2a; MIP2A; SCYB2

Mammalian Cell Puromycin

Selection:

**Vector:** pLenti-C-Myc-DDK-P2A-Puro (PS100092)

E. coli Selection: Chloramphenicol (34 ug/mL)

ORF Nucleotide The ORF insert of this clone is exactly the same as(RC204877).

Sequence:

Restriction Sites: Sgfl-Mlul

**Cloning Scheme:** 






<sup>\*</sup> The last codon before the Stop codon of the ORF.





#### Plasmid Map:



**ACCN:** NM\_002089

ORF Size: 321 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of

reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube

containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:** 1. Centrifuge at 5,000xg for 5min.

2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.

3. Close the tube and incubate for 10 minutes at room temperature.

4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid

at the bottom.

5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of

shipping when stored at -20°C.

**RefSeq:** <u>NM 002089.1</u>

RefSeq Size: 1234 bp RefSeq ORF: 324 bp



## GRO beta (CXCL2) (NM\_002089) Human Tagged Lenti ORF Clone - RC204877L3

 Locus ID:
 2920

 UniProt ID:
 P19875

 Cytogenetics:
 4q13.3

 Domains:
 IL8

**Protein Families:** Druggable Genome, Secreted Protein

Protein Pathways: Chemokine signaling pathway, Cytokine-cytokine receptor interaction, NOD-like receptor

signaling pathway

**MW:** 11.4 kDa

**Gene Summary:** This antimicrobial gene is part of a chemokine superfamily that encodes secreted proteins

involved in immunoregulatory and inflammatory processes. The superfamily is divided into four subfamilies based on the arrangement of the N-terminal cysteine residues of the mature

peptide. This chemokine, a member of the CXC subfamily, is expressed at sites of

inflammation and may suppress hematopoietic progenitor cell proliferation. [provided by

RefSeq, Sep 2014]