

Product datasheet for **RC204880L3V**

Frataxin (FXN) (NM_000144) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Frataxin (FXN) (NM_000144) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Frataxin
Synonyms:	CyaY; FA; FARR; FRDA; X25
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_000144
ORF Size:	630 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204880).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_000144.3
RefSeq Size:	7168 bp
RefSeq ORF:	633 bp
Locus ID:	2395
UniProt ID:	Q16595
Cytogenetics:	9q21.11
Domains:	Frataxin_Cyay
Protein Families:	Druggable Genome



[View online »](#)

MW: 23.1 kDa

Gene Summary: This nuclear gene encodes a mitochondrial protein which belongs to the FRATAXIN family. The protein functions in regulating mitochondrial iron transport and respiration. The expansion of intronic trinucleotide repeat GAA from 8-33 repeats to >90 repeats results in Friedreich ataxia. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2016]