

Product datasheet for **RC204881L2V**

Cathepsin G (CTSG) (NM_001911) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cathepsin G (CTSG) (NM_001911) Human Tagged ORF Clone Lentiviral Particle
Symbol:	Cathepsin G
Synonyms:	CATG; CG
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_001911
ORF Size:	765 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204881).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001911.2
RefSeq Size:	924 bp
RefSeq ORF:	768 bp
Locus ID:	1511
UniProt ID:	P08311
Cytogenetics:	14q12
Domains:	Tryp_SPC
Protein Families:	Druggable Genome, Protease



[View online »](#)

Protein Pathways: Lysosome, Neuroactive ligand-receptor interaction, Renin-angiotensin system, Systemic lupus erythematosus

MW: 28.8 kDa

Gene Summary: The protein encoded by this gene, a member of the peptidase S1 protein family, is found in azurophil granules of neutrophilic polymorphonuclear leukocytes. The encoded protease has a specificity similar to that of chymotrypsin C, and may participate in the killing and digestion of engulfed pathogens, and in connective tissue remodeling at sites of inflammation. In addition, the encoded protein is antimicrobial, with bacteriocidal activity against *S. aureus* and *N. gonorrhoeae*. Transcript variants utilizing alternative polyadenylation signals exist for this gene. [provided by RefSeq, Sep 2014]