

Product datasheet for **RC204939L3V**

LIAS (NM_006859) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	LIAS (NM_006859) Human Tagged ORF Clone Lentiviral Particle
Symbol:	LIAS
Synonyms:	HGCLAS; HUSSY-01; LAS; LIP1; LS; PDHLD
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_006859
ORF Size:	1116 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204939).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006859.2
RefSeq Size:	1764 bp
RefSeq ORF:	1119 bp
Locus ID:	11019
UniProt ID:	O43766
Cytogenetics:	4p14
Protein Pathways:	Lipoic acid metabolism, Metabolic pathways
MW:	41.9 kDa



[View online »](#)

Gene Summary:

The protein encoded by this gene belongs to the biotin and lipoic acid synthetases family. Localized in the mitochondrion, this iron-sulfur enzyme catalyzes the final step in the de novo pathway for the biosynthesis of lipoic acid, a potent antioxidant. The deficient expression of this enzyme has been linked to conditions such as diabetes, atherosclerosis and neonatal-onset epilepsy. Alternative splicing occurs at this locus, and several transcript variants encoding distinct isoforms have been identified. [provided by RefSeq, Aug 2020]