

Product datasheet for **RC204946L4V**

LRFN4 (NM_024036) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	LRFN4 (NM_024036) Human Tagged ORF Clone Lentiviral Particle
Symbol:	LRFN4
Synonyms:	FIGLER6; SALM3; SALM3.
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_024036
ORF Size:	1905 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204946).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024036.3
RefSeq Size:	2546 bp
RefSeq ORF:	1908 bp
Locus ID:	78999
UniProt ID:	Q6PJG9
Cytogenetics:	11q13.2
Domains:	LRRCT, LRR, ig, LRR_TYP, IGc2, IG, FN3
Protein Families:	Transmembrane



[View online »](#)

MW: 66.9 kDa

Gene Summary: Promotes neurite outgrowth in hippocampal neurons. May play a role in redistributing DLG4 to the cell periphery (By similarity).[UniProtKB/Swiss-Prot Function]