

Product datasheet for **RC204948L4V**

EFHD2 (NM_024329) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	EFHD2 (NM_024329) Human Tagged ORF Clone Lentiviral Particle
Symbol:	EFHD2
Synonyms:	SWS1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_024329
ORF Size:	720 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204948).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_024329.4
RefSeq Size:	2424 bp
RefSeq ORF:	723 bp
Locus ID:	79180
UniProt ID:	Q96C19
Cytogenetics:	1p36.21
Domains:	EFh
MW:	26.7 kDa



[View online »](#)

Gene Summary:

May regulate B-cell receptor (BCR)-induced immature and primary B-cell apoptosis. Plays a role as negative regulator of the canonical NF-kappa-B-activating branch. Controls spontaneous apoptosis through the regulation of BCL2L1 abundance.[UniProtKB/Swiss-Prot Function]