

Product datasheet for **RC204951A1V**

Human alpha 1 Antitrypsin (SERPINA1) (NM_001002235) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human alpha 1 Antitrypsin (SERPINA1) (NM_001002235) AAV Particle
Tag:	Myc-DDK
Symbol:	alpha 1 Antitrypsin
Synonyms:	A1A; A1AT; AAT; alpha1AT; nNIF; PI; PI1; PRO2275
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC204951 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s) TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCCGCCGCGATCGCC

ATGCCGTCTTCTGCTCTCGTGGGGCATCCTCCTGCTGGCAGGCCTGTGCTGCCTGGTCCCTGTCTCCCTGGCTGAGGATCCCCAGGGAGATGCTGCCCAGAAGACAGATACATCCCACCATGATCAGGATCACCCAACCTTCAACAAGATCACCCCAACCTGGCTGAGTTCGCCTTACGCCTATACCGCCAGCTGGCACACCAAGTCCAACAGCACAATATCTTCTTCTCCCCAGTGAGCATCGCTACAGCCTTTGCAATGCTCTCCCTGGGGACCAAGGCTGACACTCACGATGAAATCCTGGAGGGCCTGAATTTCAACCTCACGGAGATTCCGGAGGCTCAGATCCA TGAAGGCTTCCAGGAACCTCCTCCGTACCCTCAACCAGCCAGACAGCCAGCTCCAGCTGACCACCGCAATGGCCTGTTCTCAGCGAGGGCCTGAAGCTAGTGGATAAGTTTTTGGAGGATGTTAAAAAGTTGTACCACTCAGAAGCCTTCACTGTCAACTTCGGGGACACCGAAGAGGCCAAGAAACAGATCAACGATTACGTGGAGAA GGGTACTCAAGGGAAAATTGTGGATTTGGTCAAGGAGCTTGACAGAGACACAGTTTTTGTCTGCTGTAAT TACATCTTCTTTAAAGGCAAATGGGAGAGACCCTTTGAAGTCAAGGACACCGAGGAAGAGGACTTCCACG TGGACCAGGCGACCAACCGTGAAGGTGCCTATGATGAAGCGTTTAGGCATGTTTAACATCCAGCACTGTAA GAAGCTGTCCAGCTGGGTGCTGCTGATGAAATACCTGGGCAATGCCACCGCCATCTTCTCCTGCCTGAT GAGGGGAACTACAGCACCTGGAATGAAGTCAACCCAGCATATCATCAACCAAGTTCTGGAATGAAG ACAGAAGTCTGCCAGCTTACATTTACCCAACTGTCCATTACTGGAACCTATGATCTGAAGAGCGTCT GGGTCAACTGGGCATCACTAAGGTCTTCAAGCAATGGGGCTGACCTCTCCGGGTCACAGAGGAGGACCC CTGAAGCTCTCAAGGCCGTGCATAAGGCTGTGCTGACCATCGACGAGAAAGGGACTGAAGCTGTGGG CCATGTTTTTAGAGGCCATACCATGTCTATCCCCCGAGGTCAAGTTCAACAAACCTTTGTCTTCTT AATGATTGAACAAAATACCAAGTCTCCCTCTTCATGGGAAAAGTGGTGAATCCCAACCAAAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCTGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC204951 protein sequence Red=Cloning site Green=Tags(s) MPSSVSWGILLLAGLCLVPVSLAEDPQGDAQKTDTSHHDDHPTFNKITPNLAFAFSLYRQLAHQSN STNIFFSPVSIATAFAMLSLGTKADTHDEILEGLNFNLTEIPEAQIHEGFQELLRTLNPDSQLQLTTGN GLFLSEGLKLVDFLEDVKKLYHSEFTVNFGDTEEAKKQINDYVEKGTQGKIVDLVKELDRDVFALVN YIFFKKGWERPFVEVKDTEEDFHVDQATTVKVPMKRLGMFNIQHCKLSSWVLLMKYLGNATAIFFLPD EGKLQHLENELTHDIITKFLNEDRRSASLHLPKLSITGTYDLKSVLGQLGITKVFSGADLSGVTEEAP LKLSKAVHKAVLTIDEKGTEAAGAMFLEAIPMSIPPEVKFNKPFVFLMIEQNTKSPLFMGKVVNPTQK TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_001002235
ORF Size:	1254 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001002235.1
RefSeq Size:	3199 bp
RefSeq ORF:	1257 bp
Locus ID:	5265
UniProt ID:	P01009
Cytogenetics:	14q32.13
MW:	46.7 kDa