

## Product datasheet for **RC204995L3V**

### ETV3 (NM\_005240) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	ETV3 (NM_005240) Human Tagged ORF Clone Lentiviral Particle
Symbol:	ETV3
Synonyms:	METS; PE-1; PE1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_005240
ORF Size:	429 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204995).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_005240.1</a>
RefSeq Size:	1443 bp
RefSeq ORF:	432 bp
Locus ID:	2117
UniProt ID:	<a href="#">P41162</a>
Cytogenetics:	1q23.1
Protein Families:	Transcription Factors
MW:	16.9 kDa



[View online »](#)

**Gene Summary:**

Transcriptional repressor that contribute to growth arrest during terminal macrophage differentiation by repressing target genes involved in Ras-dependent proliferation. Represses MMP1 promoter activity.[UniProtKB/Swiss-Prot Function]