

Product datasheet for **RC204997L4V**

C9orf78 (NM_016520) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	C9orf78 (NM_016520) Human Tagged ORF Clone Lentiviral Particle
Symbol:	C9orf78
Synonyms:	bA409K20.3; CSU2; HCA59; HSPC220
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_016520
ORF Size:	867 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC204997).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_016520.2
RefSeq Size:	1835 bp
RefSeq ORF:	870 bp
Locus ID:	51759
UniProt ID:	Q9NZ63
Cytogenetics:	9q34.11
MW:	33.7 kDa



[View online »](#)

Gene Summary: Involved in the regulation of telomeric heterochromatin assembly and control of telomere length.[UniProtKB/Swiss-Prot Function]