

## Product datasheet for **RC205002L3V**

### SLIRP (NM\_031210) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	SLIRP (NM_031210) Human Tagged ORF Clone Lentiviral Particle
Symbol:	SLIRP
Synonyms:	C14orf156; DC50; PD04872
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_031210
ORF Size:	327 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205002).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_031210.3</a> , <a href="#">NP_112487.1</a>
RefSeq Size:	420 bp
RefSeq ORF:	330 bp
Locus ID:	81892
UniProt ID:	<a href="#">Q9GZT3</a>
Cytogenetics:	14q24.3
Domains:	RRM
MW:	12.3 kDa



[View online »](#)

**Gene Summary:**

Steroid receptor RNA activator (SRA, or SRA1; MIM 603819) is a complex RNA molecule containing multiple stable stem-loop structures that functions in coactivation of nuclear receptors. SLIRP interacts with stem-loop structure-7 of SRA (STR7) and modulates nuclear receptor transactivation (Hatchell et al., 2006 [PubMed 16762838]).[supplied by OMIM, Mar 2008]