

Product datasheet for **RC205010L4V**

SDHC (NM_003001) Human Tagged ORF Clone Lentiviral Particle

Product data:

| | |
|---------------------------|--|
| Product Type: | Lentiviral Particles |
| Product Name: | SDHC (NM_003001) Human Tagged ORF Clone Lentiviral Particle |
| Symbol: | SDHC |
| Synonyms: | CYB560; CYBL; PGL3; QPS1; SDH3 |
| Mammalian Cell Selection: | Puromycin |
| Vector: | pLenti-C-mGFP-P2A-Puro (PS100093) |
| Tag: | mGFP |
| ACCN: | NM_003001 |
| ORF Size: | 507 bp |
| ORF Nucleotide Sequence: | The ORF insert of this clone is exactly the same as(RC205010). |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| RefSeq: | NM_003001.2 |
| RefSeq Size: | 2858 bp |
| RefSeq ORF: | 510 bp |
| Locus ID: | 6391 |
| UniProt ID: | Q99643 |
| Cytogenetics: | 1q23.3 |
| Domains: | Sdh_cyt |
| Protein Families: | Druggable Genome, Transmembrane |




[View online »](#)

Protein Pathways: Alzheimer's disease, Citrate cycle (TCA cycle), Huntington's disease, Metabolic pathways, Oxidative phosphorylation, Parkinson's disease

MW: 18.6 kDa

Gene Summary: This gene encodes one of four nuclear-encoded subunits that comprise succinate dehydrogenase, also known as mitochondrial complex II, a key enzyme complex of the tricarboxylic acid cycle and aerobic respiratory chains of mitochondria. The encoded protein is one of two integral membrane proteins that anchor other subunits of the complex, which form the catalytic core, to the inner mitochondrial membrane. There are several related pseudogenes for this gene on different chromosomes. Mutations in this gene have been associated with paragangliomas. Alternatively spliced transcript variants have been described. [provided by RefSeq, May 2013]

Product images:



[RC205010L4] was used to prepare Lentiviral particles using [TR30037] packaging kit. HEK293T cells were transduced with RC205010L4V particle to overexpress human SDHC-mGFP fusion protein.