

## Product datasheet for **RC205025L2V**

### TIMM17A (NM\_006335) Human Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	TIMM17A (NM_006335) Human Tagged ORF Clone Lentiviral Particle
Symbol:	TIMM17A
Synonyms:	TIM17; TIM17A
Mammalian Cell Selection:	None
Vector:	pLenti-C-mGFP (PS100071)
Tag:	mGFP
ACCN:	NM_006335
ORF Size:	513 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205025).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_006335.1</a>
RefSeq Size:	1675 bp
RefSeq ORF:	516 bp
Locus ID:	10440
UniProt ID:	<a href="#">Q99595</a>
Cytogenetics:	1q32.1
Domains:	Tim17
Protein Families:	Transmembrane



[View online »](#)

**MW:** 18 kDa

**Gene Summary:** Essential component of the TIM23 complex, a complex that mediates the translocation of transit peptide-containing proteins across the mitochondrial inner membrane.  
[UniProtKB/Swiss-Prot Function]