

Product datasheet for **RC205039**

PIH1D2 (NM_138789) Human Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	PIH1D2 (NM_138789) Human Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	PIH1D2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>RC205039 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGACATCCTCAAAAGGTCTGCTTACCCAAGTTACTCAGTTTTGGAACCTCCTAGATGATCTAGCTC
AGAGTGACCCTGAGGGCTATGAGAAGTTTATTCAGCAGCAGCTGAAAGAAGGGAAACAGCTCTGTGCTGC
CCCAGAACCACAGCTTTGTCTACAGACCAGGATTCTGAAACCAAAAGAAAAATACTTTTATCAACCTG
TGTCAGTGGACAAGGATCCCAGCTCCCAATCAACCACTCATCCAGTACCTCTAACTGTTGGCAAACAG
AAGATACAACTGAGATATCAGATGCTTACACAGTCATTGATGTTGCCTACAATCCTGATGTTCTTCATGC
AGCAGAAAAGGACCAAGTGAAGAAAAATCAGTTAATTCAGATGGCCATGAAATGCATTGAGGAGAAATTC
CAGTTCACCTCTCACACTCTTACCATATTACCAATTTAGAATAAAAGGAAGCATTCAAAGAATGAAAC
AAAATCTGATGGGAATCCAACTGATTCCATAGATTTAAGAGAAAAATGAGAAGGGAACAACTCTTGG
ACAGATACGAAGCAGTACTATGAGCAATCCAGATCACTTTCCTCAACTGTTACTGCCAAAAGACCAAGTT
TCAGGCAAAGCAGTGTGCTGATAGAAGAGATTCCAGTACTGAAATCCAGGTGGAGATGAAGATGCCAG
CCTATGAACTAAAAATTGTGCATGATCACAGTGAGAAACCTCTGAAAATTGAGTTGAAAGTTGAATTACC
TGGTATTAATTCTGTCTCTCTGTGACCTTAGTGTCTGAGGATGATTTATTGATTGAAGTATCTGAG
AAGTACAGATTACATCTGAATCTTCCAAAACCTATTGATACTGAAATGACCACAGCAAAATTTATCAAAG
AAAAATCCACGCTAATCATCACAATGCCTTTGGTG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA


[View online »](#)

Protein Sequence: >RC205039 protein sequence
 Red=Cloning site Green=Tags(s)


METSSKGLLTQVTQFWNLLDDLAQSDPEGYEKFIQQQLKEGKQLCAAPEPQLCLQTRILKPKEKILFINL
 CQWTRIPAPQSTTHPVPLTVGKPEDTTEISDAYTVIDVAYNPDLHAAEKDQVKKNQLIQMAMKCIEEF
 QFTLSHSYHITKFRIKGSIQRMKQNLMGITQDSIDLREKMRRELTLGQIRSSTMSNPDHFPQLLLPKDQV
 SGKAVCLIEEISSTEIQVEMKMPAYELKIVHDHSEKPLKIELKVELPGINSVSLCDLSVSEDDLLIEVSE
 KYRLHLNLPKLIDTEMTAKFIKEKSTLIITMPLV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6421_e06.zip

Restriction Sites: SgfI-MluI

Cloning Scheme:



ACCN: NM_138789

ORF Size: 945 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_138789.3](#)

RefSeq Size: 1252 bp

RefSeq ORF: 948 bp


Locus ID: 120379

UniProt ID: [Q8WWB5](#)


Cytogenetics: 11q23.1

MW: 36 kDa


Product images:




Circular map for RC205039



HEK293T cells were transfected with the pCMV6-ENTRY control (Cat# [PS100001], Left lane) or pCMV6-ENTRY PIH1D2 (Cat# RC205039, Right lane) cDNA for 48 hrs and lysed. Equivalent amounts of cell lysates (5 ug per lane) were separated by SDS-PAGE and immunoblotted with anti-PIH1D2 (Cat# [TA502467]). Positive lysates [LY408526] (100ug) and [LC408526] (20ug) can be purchased separately from OriGene.



Western blot validation of overexpression lysate (Cat# [LY408526]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC205039 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified PIH1D2 protein (Cat# [TP305039]). The protein was produced from HEK293T cells transfected with PIH1D2 cDNA clone (Cat# RC205039) using MegaTran 2.0 (Cat# [TT210002]).