

Product datasheet for **RC205042L1V**

IFI44 (NM_006417) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	IFI44 (NM_006417) Human Tagged ORF Clone Lentiviral Particle
Symbol:	IFI44
Synonyms:	MTAP44; p44; TLDC5
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_006417
ORF Size:	1332 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205042).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_006417.2
RefSeq Size:	1742 bp
RefSeq ORF:	1335 bp
Locus ID:	10561
UniProt ID:	Q8TCB0
Cytogenetics:	1p31.1
MW:	50.5 kDa
Gene Summary:	This protein aggregates to form microtubular structures.[UniProtKB/Swiss-Prot Function]



[View online »](#)