

Product datasheet for **RC205064**

H4-16 (NM_175054) Human Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: H4-16 (NM_175054) Human Tagged ORF Clone
Tag: Myc-DDK
Symbol: H4-16
Synonyms: H4/p; H4C1; H4C2; H4C3; H4C4; H4C5; H4C6; H4C8; H4C9; H4C11; H4C12; H4C13; H4C14; H4C15; HIST4H4
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >RC205064 representing NM_175054
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTCTGGGCGAGGTAAAGGTGGCAAGGGGCTGGGTAAGGGAGGCGCAAGCGCCACCGGAAGGTGCTGC
GGGACAATATCCAAGGCATTACAAAGCCGGCGATTGCGCGTCTCGCCCGACGTGGGGCGTCAAGCGCAT
TTCTGGTCTCATCTACGAGGAGACCCGGGGAGTCTCAAAGTCTTCTGGAGAACGTGATCCGTGACGCG
GTGACTTACACGGAGCACGCCAAGCGCAAGACCGTCACGGCCATGGATGTGGTGTACGCGCTGAAACGCC
AGGGTCGCACCCTTTATGGTTTCGGCGGT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RC205064 representing NM_175054
Red=Cloning site Green=Tags(s)
MSGRGKGGKGLGKGGAKRHRKVLRDNIQGITKPAIRRLARRGGVKRISGLIYEETRGVLKVFLENVIRDA
VTYTEHAKRKTVTAMDVVYALKRQGRPLYGFGG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Chromatograms: https://cdn.origene.com/chromatograms/mk6439_f08.zip

Restriction Sites: Sgfl-Mlul



[View online »](#)

Cloning Scheme:


ACCN: NM_175054

ORF Size: 309 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: [NM_175054.2](#), [NP_778224.1](#)

RefSeq Size: 412 bp


RefSeq ORF: 312 bp

Locus ID: 121504
UniProt ID: [P62805](#)
Cytogenetics: 12p12.3


Protein Pathways: Systemic lupus erythematosus
MW: 11.2 kDa

Gene Summary: Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. Nucleosomes consist of approximately 146 bp of DNA wrapped around a histone octamer composed of pairs of each of the four core histones (H2A, H2B, H3, and H4). The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the DNA between the nucleosomes to form higher order chromatin structures. This gene is intronless and encodes a replication-dependent histone that is a member of the histone H4 family. Transcripts from this gene lack polyA tails; instead, they contain a palindromic termination element. [provided by RefSeq, Aug 2015]


Product images:



Circular map for RC205064



Western blot validation of overexpression lysate (Cat# [LY406378]) using anti-DDK antibody (Cat# [TA50011-100]). Left: Cell lysates from untransfected HEK293T cells; Right: Cell lysates from HEK293T cells transfected with RC205064 using transfection reagent MegaTran 2.0 (Cat# [TT210002]).



Coomassie blue staining of purified HIST4H4 protein (Cat# [TP305064]). The protein was produced from HEK293T cells transfected with HIST4H4 cDNA clone (Cat# RC205064) using MegaTran 2.0 (Cat# [TT210002]).