

Product datasheet for **RC205086A1V**

Human Parvin alpha (PARVA) (NM_018222) AAV Particle

Product data:

Product Type:	AAV Particles
Product Name:	Human Parvin alpha (PARVA) (NM_018222) AAV Particle
Tag:	Myc-DDK
Symbol:	Parvin alpha
Synonyms:	CH-ILKBP; MXRA2
Mammalian Cell Selection:	None
Vector:	pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence:	>RC205086 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGC**

ATGGCCACCTCCCCGAGAAGTCGCCTTCTGTCCCAAGTCTCCCACTCCCAAGTCGCCCCGTCCCGCA
AGAAAGATGATTCCTTCTGGGAAACTCGGAGGGACCCTGGCCCGGAGGAAGAAAGCCAAGGAGGTGTC
CGAGCTGCAGGAGGAGGAATGAACGCCATCAACCTGCCCTCAGCCCAATTCCCTTTGAGCTGGACCCC
GAGGACACGATGCTGGAGGAGAATGAGGTGCGAACAATGGTGGATCCAACTCACGCAGTGACCCCAAGC
TTCAAGAAGTGAAGGTATTAATTGACTGGATTAATGATGTGTTGGTTGGAGAAAGAATCATTGTGAA
AGACCTAGCTGAAGATTTGTATGATGGACAAGTCCTGCAGAAGCTTTTCGAGAAACTGGAGAGTGAGAAG
CTAAATGTGGCTGAGGTACCCAGTCAGAGATTGCTCAGAAGCAAAACTGCAGACTGTCTGGAGAAGA
TCAATGAAACCTGAACTTCCTCCCAGGAGCATCAAGTGGAATGTGGATTCTGTTATGCCAAGAGCCT
GGTGGCCATCTTACACCTGCTCGTTGCTCTGTCTCAGTATTTCCGCGCACCAATTCGACTCCCAGACCAT
GTTTCCATCCAAGTGTTGTGGTCCAGAAACGAGAAGGAATCCTCCAGTCTCGGCAAATCCAAGAGGAAA
TAACTGGTAACACAGAGGCTCTTCCGGGAGGCATGAACGTGATGCCTTTGACACCTTGTTGCACCATGC
CCCAGACAAGCTGAATGTGGTAAAAAGACACTCATCACTTTCTGTAACAAGCACCTGAATAAACTGAAC
CTGGAGGTCACAGAAGTGGAAACCCAGTTGCAGATGGGGTGTACCTGGTGTCTCATGGGGCTCCTGG
AGGGCTACTTTGTGCCCTGCACAGCTTCTTCCCTGACCCCGGACAGCTTTGAACAGAAGGTCTTGAATGT
CTCCTTTGCCTTTGAGCTCATGCAAGATGGAGGTTGGAAAAGCCAAACCGCGGCCAGAAGACATAGTC
AACTGTGACCTGAAATCTACACTACGAGTGTTGTACAACCTCTTCACCAAGTACCGTAACGTGGAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>RC205086 protein sequence Red=Cloning site Green=Tags(s) MATSPQKSPSPKSPTPKSPPSRKDDSFGLGKGGTLARRKKAKEVSELQEEGMNAINLPLSPIPFELDP EDTMLEENEVRTMVDPNRSRDPKLQELMKVLIDWINDVLVGERIIVKDLAEDLYDGQVLQKLFKLESEK LNVAEVTQSEIAQKQKLQTVLEKINETLKLPPRSIKWNVDSVHAKSLVAILHLLVALSQYFRAPIRLPDH VSIQVVVVQKREGILQSRQIQEEITGNTEALSGRHERDAFDLFDHAPDKLNVVKKTLITFVNKHLNKLN LEVTELETQFADGVYLVLLMGLLEGYFVPLHSFFLTPDSFEQKVLNVSAFELMQDGGLEKPKPRPEDIV NCDLKSTLRVLYNLFTKYRNV TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Species:	Human
Serotype:	AAV-2
ACCN:	NM_018222
ORF Size:	1116 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_018222.2
RefSeq Size:	8719 bp
RefSeq ORF:	1119 bp
Locus ID:	55742
UniProt ID:	Q9NVD7
Cytogenetics:	11p15.3
MW:	42.2 kDa