

## Product datasheet for **RC205093L4V**

### **MVB12B (NM\_001011703) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	MVB12B (NM_001011703) Human Tagged ORF Clone Lentiviral Particle
Symbol:	MVB12B
Synonyms:	C9orf28; FAM125B
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001011703
ORF Size:	663 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205093).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001011703.1</a>
RefSeq Size:	2705 bp
RefSeq ORF:	666 bp
Locus ID:	89853
UniProt ID:	<a href="#">Q9H7P6</a>
Cytogenetics:	9q33.3
Protein Pathways:	Endocytosis
MW:	24.5 kDa



[View online »](#)

**Gene Summary:**

The protein encoded by this gene is a component of the ESCRT-I complex, a heterotetramer, which mediates the sorting of ubiquitinated cargo protein from the plasma membrane to the endosomal vesicle. ESCRT-I complex plays an essential role in HIV budding and endosomal protein sorting. Depletion and overexpression of this and related protein (MVB12A) inhibit HIV-1 infectivity and induce unusual viral assembly defects, indicating a role for MVB12 subunits in regulating ESCRT-mediated virus budding. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2011]