

Product datasheet for **RC205104L4V**

C1orf150 (GCSAML) (NM_145278) Human Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	C1orf150 (GCSAML) (NM_145278) Human Tagged ORF Clone Lentiviral Particle
Symbol:	C1orf150
Synonyms:	C1orf150
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_145278
ORF Size:	405 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205104).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_145278.1
RefSeq Size:	3863 bp
RefSeq ORF:	408 bp
Locus ID:	148823
UniProt ID:	Q5JQS6
Cytogenetics:	1q44
MW:	15.7 kDa



[View online »](#)

Gene Summary:

This gene encodes a protein thought to be a signaling molecule associated with germinal centers, the sites of proliferation and differentiation of mature B lymphocytes. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Aug 2013]