

## Product datasheet for **RC205111L4V**

### **RMI2 (NM\_152308) Human Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	RMI2 (NM_152308) Human Tagged ORF Clone Lentiviral Particle
Symbol:	RMI2
Synonyms:	BLAP18; C16orf75
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_152308
ORF Size:	441 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RC205111).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_152308.1</a>
RefSeq Size:	1462 bp
RefSeq ORF:	444 bp
Locus ID:	116028
UniProt ID:	<a href="#">Q96E14</a>
Cytogenetics:	16p13.13
MW:	15.7 kDa



[View online »](#)

**Gene Summary:**

RM12 is a component of the BLM (RECQL3; MIM 604610) complex, which plays a role in homologous recombination-dependent DNA repair and is essential for genome stability (Xu et al., 2008 [PubMed 18923082]).[supplied by OMIM, Nov 2008]